

A-DHS ANTIROTATOR DYNAMIC HIP SCREW PLATES

Indications: Subtrochanteric, pertrochanteric, intertrochanteric femoral fractures and collum femoris fractures are suitable for applying A-DHS.
Antirotator Dynamic Hip Screw and Plate System, which has been improved by performing significant revisions to DHS System, strengthens rotational stability of femoral head owing to the unique blade design. The blade can be sent over the lag screw and also allows sliding compression. A-DHS consists of four basic parts; Plate, Lag Screw, Antirotator Blade and Compression Screw.

Plate Profile: Thickness: 8.0 mm, Width: 19.0 mm, Distance between holes: 16.0 mm

Cat. Code	Barcode	Description	Material	Length (mm)
33227135003	8699931025378	A-DHS LOCK PLATE 135° TI 3 HOLES	Ti	73
33227135004	8699931025385	A-DHS LOCK PLATE 135° TI 4 HOLES	Ti	89
33227135005	8699931025392	A-DHS LOCK PLATE 135° TI 5 HOLES	Ti	105
33227135006	8699931025408	A-DHS LOCK PLATE 135° TI 6 HOLES	Ti	127
33227135007	8699931025415	A-DHS LOCK PLATE 135° TI 7 HOLES	Ti	143
33227135008	8699931025422	A-DHS LOCK PLATE 135° TI 8 HOLES	Ti	159
33227135009	8699931025439	A-DHS LOCK PLATE 135° TI 9 HOLES	Ti	175
33227135010	8699931025309	A-DHS LOCK PLATE 135° TI 10 HOLES	Ti	191
33227135011	8699931025316	A-DHS LOCK PLATE 135° TI 11 HOLES	Ti	207



A-DHS LAG SCREW -Titanium			
Cat. Code	Barcode	Length	Qty
20720600001	8699931025231	60 MM	1
20720650001	8699931025248	65 MM	1
20720700001	8699931025255	70 MM	1
20720750001	8699931025590	75 MM	1
20720800001	8699931025262	80 MM	1
20720850001	8699931025279	85 MM	1
20720900001	8699931025286	90 MM	1
20720950001	8699931025293	95 MM	1
20720100001	8699931025163	100 MM	1
20720105001	8699931025170	105 MM	1
20720110001	8699931025187	110 MM	1



ANTIROT. BLADE (A-DHS) -Titanium			
Cat. Code	Barcode	Length	Qty
21620000060	8699931025514	60 MM	1
21620000065	8699931025521	65 MM	1
21620000070	8699931025538	70 MM	1
21620000075	8699931025545	75 MM	1
21620000080	8699931025552	80 MM	1
21620000085	8699931025569	85 MM	1
21620000090	8699931025576	90 MM	1
21620000095	8699931025583	95 MM	1
21620000100	8699931025446	100 MM	1
21620000105	8699931025453	105 MM	1
21620000110	8699931025460	110 MM	1



A-DHS COMPRESSION SCREW -Titanium			
Cat. Code	Barcode	Length	Qty
21027000012	8699931025804	25 mm	11



Recommended Screws

- 4.5 mm Cortical Screw
- 4.5 mm Lock Cortical Screw
- 6.5 mm Cancellous Screw
- 6.5 mm Lock Cancellous Screw



LW-06 A-DHS PLATES SET (Screws)

Row	Cat. Code	Barcode	Description	Quantity
1	0201200	8699931010817	SCREW TRAY	1
2	20227160045	8698673417793	CORT.LOCK.SCREW Ø4.5 MM (16mm-50mm 2mm incr) 5 ea.	90
3	20220160045	8698673417298	CORT. SCREW Ø4.5 MM (16mm-50mm 2mm incr) 4 ea.	72
4	20427300065	8698673416956	CANCEL.LOCK.SCREW Ø6.5 MM (30mm-100mm 5mm incr) 30-35mm 3 ea.	6
	20427400065	8698673416970	CANCEL.LOCK.SCREW Ø6.5 MM 40-45-50-60mm 6 ea.	24
	20427550065	8698673417007	CANCEL.LOCK.SCREW Ø6.5 MM 55-65-70-75-80-85-90-95-100mm 3ea.	27
5	20420303265	8698673421653	CANCELLOUS SCREW (32 THREAD) 6.5 MM (30mm-100mm 2mm incr) 3 ea.	30


Recommended Screws

Ø 4.5 mm.Locking Cortical

Ø 4.5 mm Cortical

Ø 6.5 mm Locking Cortical

Ø 6.5 mm Cancellous



LOCKING CORTICAL SCREW Ø 4.5 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20227160045	8698673417793	20217160045	8698673449305	16	5
20227180045	8698673417809	20217180045	8698673449312	18	5
20227200045	8698673417816	20217200045	8698673449329	20	5
20227220045	8698673417823	20217220045	8698673449336	22	5
20227240045	8698673417830	20217240045	8698673449343	24	5
20227260045	8698673417847	20217260045	8698673449350	26	5
20227280045	8698673417854	20217280045	8698673449367	28	5
20227300045	8698673417861	20217300045	8698673449374	30	5
20227320045	8698673417878	20217320045	8698673449381	32	5
20227340045	8698673417885	20217340045	8698673449398	34	5
20227360045	8698673417892	20217360045	8698673449404	36	5
20227380045	8698673417908	20217380045	8698673449411	38	5
20227400045	8698673417915	20217400045	8698673449428	40	5
20227420045	8698673417922	20217420045	8698673449435	42	5
20227440045	8698673417939	20217440045	8698673449442	44	5
20227460045	8698673417946	20217460045	8698673449459	46	5
20227480045	8698673417953	20217480045	8698673449466	48	5
20227500045	8698673417960	20217500045	8698673449473	50	5

CORTICAL SCREW Ø 4.5 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20220160045	8698673417298	20210160045	8698673412538	16	4
20220180045	8698673417304	20210180045	8698673412545	18	4
20220200045	8698673417311	20210200045	8698673412552	20	4
20220220045	8698673417328	20210220045	8698673412569	22	4
20220240045	8698673417335	20210240045	8698673412576	24	4
20220260045	8698673417342	20210260045	8698673412583	26	4
20220280045	8698673417359	20210280045	8698673412590	28	4
20220300045	8698673417366	20210300045	8698673412606	30	4
20220320045	8698673417373	20210320045	8698673412613	32	4
20220340045	8698673417380	20210340045	8698673412620	34	4
20220360045	8698673417397	20210360045	8698673412637	36	4
20220380045	8698673417403	20210380045	8698673412644	38	4
20220400045	8698673417410	20210400045	8698673412651	40	4
20220420045	8698673417427	20210420045	8698673412668	42	4
20220440045	8698673417434	20210440045	8698673412675	44	4
20220460045	8698673417441	20210460045	8698673412682	46	4
20220480045	8698673417458	20210480045	8698673412699	48	4
20220500045	8698673417465	20210500045	8698673412705	50	4

LOCKING CANCELLOUS SCREW Ø 6.5 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20427300065	8698673416956	20417300065	8698673449749	30	3
20427350065	8698673416963	20417350065	8698673449756	35	3
20427400065	8698673416970	20417400065	8698673449763	40	6
20427450065	8698673416987	20417450065	8698673449770	45	6
20427500065	8698673416994	20417500065	8698673459793	50	6
20427550065	8698673417007	20417550065	8698673449602	55	3
20427600065	8698673417014	20417600065	8698673449619	60	6
20427650065	8698673417021	20417650065	8698673449626	65	3
20427700065	8698673417038	20417700065	8698673449633	70	3
20427750065	8698673417045	20417750065	8698673449640	75	3
20427800065	8698673417052	20417800065	8698673449657	80	3
20427850065	8698673422391	20417850065	8698673449664	85	3
20427900065	8698673424654	20417900065	8698673449671	90	3
20427950065	8698673424661	20417950065	8698673449688	95	3
20427100065	8698673476691	20417100065	8698673449695	100	3

CANCELLOUS SCREW (32 THREAD) Ø 6.5 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20420303265	8698673421653	20410303265	8698673411579	30	2
20420353265	8698673421660	20410353265	8698673411586	35	2
20420403265	8698673421677	20410403265	8698673411593	40	2
20420453265	8698673421684	20410453265	8698673411609	45	2
20420503265	8698673421691	20410503265	8698673411616	50	2
20420553265	8698673421707	20410553265	8698673411623	55	2
20420603265	8698673421714	20410603265	8698673411630	60	2
20420653265	8698673421721	20410653265	8698673411647	65	2
20420703265	8698673421738	20410703265	8698673411654	70	2
20420753265	8698673421745	20410753265	8698673411661	75	2
20420803265	8698673421752	20410803265	8698673411678	80	2
20420853265	8698673421769	20410853265	8698673411685	85	2
20420903265	8698673421776	20410903265	8698673411692	90	2
20420953265	8698673421783	20410953265	8698673411708	95	2
20421003265	8698673421622	20411003265	8698673411531	100	2

LW-06 A DHS PLATES SET (Implants & Instruments)

Row	Cat. Code	Barcode	Description	Quantity
1	0207100	8699931017212	1. DESIGN TRAY	1
2	33227135003	8699931025378	PLATES A DHS 135° (3,4,5,6,7,8,9,10,11,12 holes)	9
3	03325100000	8698673495791	DHS ADJUSTABLE ANGLE GUIDE	1
4	01193000023	8698673493780	T QUICK HANDLE	1
5	0207200	8699931017229	2. DESIGN TRAY	1
6	20720600001	8699931025231	A-DHS LAG SCREW Ti (60 mm-110 mm 5 mm incr.) ea.1	11
7	21620000060	8699931025514	ANTIROT. BLADE (A-DHS) (60 mm-110 mm 5 mm incr.) ea 1	11
8	21027000012	8699931025804	A DHS COMPRESSION SCREW Ti	11
9	01017013013	8699931026061	COUPLING SCREW FOR A DHS	1
10	03325600001	8699931026092	T HANDLE WRENCH A-DHS	1
11	22310160035	8698673458826	BONE DRILL 3.5X60X160 MM	2
12	01210030045	8699931026672	GRADUATED DRILL BIT Ø4.5 MMX300 MM (STOPPER)	1
13	04551208260	8699931029628	K-WIRE TUBE Ø12XØ8X260 MM	1
14	23412250025	8698673452770	KIRSCHNER WIRE TREADED POINT 2.5X250 MM	2
15	21510250140	8698673454095	STEINMANN PINS 4X250 MM	1
16	02171000100	8699931025996	A-DHS ANTIROT.BLADE INSERTER	1
17	01017061015	8699931026016	A-DHS BLADE GUIDE	1
18	01017018014	8699931026009	A-DHS BLADE DRILL SLEEVE	1
19	02050009020	8680858410989	LENGTH MEASURING DEVICE DHS & A-DHS	1
20	0207300	8699931017236	3. DESIGN TRAY	1
21	02070000001	8698673440883	HAMMER (MISS A CHS)	1
22	03325800000	8698673495357	GUIDE SHAFT	1
23	02005000350	8698673490765	DRILL SLEEVE GUIDE Ø 3.5 MM	2
24	02005000450	8698673495814	DRILL SLEEVE GUIDE 4.5 MM	2
25	03325500000	8698673493773	TAP FOR LAG SCREW	1
26	08040000151	8698673494800	PLATE IMPACTOR	1
27	02005000100	8698673497047	DEPTH GAUGE - LOWER EXTREMITY PLATES	1
28	02010101002	8698673493308	SOFT SCREW DRIVER QUICK LARGE	1
29	02020101035	8698673493322	SCREW DRIVER QUICK Ø 3.5 MM	1
30	02030101045	8698673440715	QUICK TAP Ø 4.5 MM	1
31	03325400000	8698673493568	TRIPLE REAMER FOR DHS	1
	00560270150	8699931010275	CONTAINER 560X270X150 MM	1

