

tst | Orthopedic Implants

LW-01 SS
LW-02 Ti
CABLE
PLATES
SET

LW-01 SS / LW-02 Ti Cable Plates Sets


CABLE PLATES

Indications: Hip and acetabular fractures, Prophylactic banding in total joint replacements and reattachment of the greater trochanter following osteotomy in total hip arthroplasty or fractures.

Stainless steel plates can be used with stainless steel (CrNi) cables. Titanium plates can be used with both CrNi and CoCrMo cables.

	TST Code	Barcode	Description	Side	Material	Length (mm)
LW-01 SS	36810003001	8699931027730	PROX. FEMUR LOCK. HOOK PLATE SMALL	-	S.STEEL	45
	36810002002	8699931027747	PROX. FEMUR LOCK. HOOK PLATE MEDIUM	-	S.STEEL	60
	36810000130	8699931027754	PROX. FEMUR LOCK. HOOK PLATE 3 HOLES	-	S.STEEL	130
	36815000180	8699931027761	PROX. FEMUR LOCK. HOOK PLATE 5 HOLES	R	S.STEEL	180
	36815000255	8699931027778	PROX. FEMUR LOCK. HOOK PLATE 8 HOLES	R	S.STEEL	255
	36816000180	8699931027792	PROX. FEMUR LOCK. HOOK PLATE 5 HOLES	L	S.STEEL	180
	36816000255	8699931027808	PROX. FEMUR LOCK. HOOK PLATE 8 HOLES	L	S.STEEL	255
	36817000040	8699931028218	DIAPHYSEAL LOCK. CABLE PLATE 4 HOLES	-	S.STEEL	137
	36817000060	8699931028225	DIAPHYSEAL LOCK. CABLE PLATE 6 HOLES	-	S.STEEL	187
	36817000080	8699931028232	DIAPHYSEAL LOCK. CABLE PLATE 8 HOLES	-	S.STEEL	237
	36817000100	8699931028249	DIAPHYSEAL LOCK. CABLE PLATE 10 HOLES	-	S.STEEL	287
	26810000075	8699931027839	CABLE CLAMP SCREW		S.STEEL	
	36810000004	8699931027822	CABLE CLAMP		S.STEEL	
	30110000002	8680858420636	CABLE+SWAGE		CrNi	760
LW-02 Ti	36820003001	8699931027655	PROX. FEMUR LOCK. HOOK PLATE SMALL	-	Ti	45
	36820002002	8699931027662	PROX. FEMUR LOCK. HOOK PLATE MEDIUM	-	Ti	60
	36820000130	8699931027679	PROX. FEMUR LOCK. HOOK PLATE 3 HOLES	-	Ti	130
	36825000180	8699931027686	PROX. FEMUR LOCK. HOOK PLATE 5 HOLES	R	Ti	180
	36825000255	8699931027693	PROX. FEMUR LOCK. HOOK PLATE 8 HOLES	R	Ti	255
	36826000180	8699931027716	PROX. FEMUR LOCK. HOOK PLATE 5 HOLES	L	Ti	180
	36826000255	8699931027723	PROX. FEMUR LOCK. HOOK PLATE 8 HOLES	L	Ti	255
	36827000040	8699931028256	DIAPHYSEAL LOCK. CABLE PLATE 4 HOLES	-	Ti	137
	36827000060	8699931028263	DIAPHYSEAL LOCK. CABLE PLATE 6 HOLES	-	Ti	187
	36827000080	8699931028270	DIAPHYSEAL LOCK. CABLE PLATE 8 HOLES	-	Ti	237
	36827000100	8699931028287	DIAPHYSEAL LOCK. CABLE PLATE 10 HOLES	-	Ti	287
	26820000075	8699931028102	CABLE CLAMP SCREW		Ti	
	36820000004	8699931028058	CABLE CLAMP		Ti	
	30110000002	8680858420636	CABLE+SWAGE		CrNi	760
30110000004	8680858420643	CABLE+SWAGE		CoCrMo	760*	


*opsiyonel



Recommended Screws for LW-01 5.0 mm Cortical S. Steel, 5.0 mm Lock. Cortical S. Steel, 5.5 mm Cancellous S. Steel, 5.5 mm Lock. Cancellous S. Steel
Recommended Screws for LW-02 5.0 mm Cortical, 5.0 mm Lock. Cortical, 5.5 mm Cancellous, 5.5 mm Lock. Cancellous


LW-01 SS / LW-02 Ti Cable Plates Sets

Screws




LOCK CORTEX SCREW CONICAL HEAD SELFT Ø 5 MM

Stainless Steel		Titanium		Length	Qty
TST Code	Barcode	TST Code	Barcode		
20117120050	8698673461567	20127120050	8698673476394	12	3
20117140050	8698673461574	20127140050	8698673445185	14	3
20117160050	8698673446472	20127160050	8698673429598	16	3
20117180050	8698673446489	20127180050	8698673429604	18	3
20117200050	8698673446496	20127200050	8698673429611	20	3
20117220050	8698673446502	20127220050	8698673429628	22	3
20117240050	8698673446519	20127240050	8698673429635	24	3
20117260050	8698673446526	20127260050	8698673429642	26	3
20117280050	8698673446533	20127280050	8698673429659	28	3
20117300050	8698673446540	20127300050	8698673429666	30	3
20117320050	8698673446557	20127320050	8698673429673	32	3
20117340050	8698673446564	20127340050	8698673429680	34	3
20117360050	8698673446571	20127360050	8698673429697	36	3
20117380050	8698673446588	20127380050	8698673429703	38	3
20117400050	8698673446595	20127400050	8698673429710	40	3
20117420050	8698673461581	20127420050	8698673445192	42	3
20117440050	8698673461598	20127440050	8698673445208	44	3
20117460050	8698673461604	20127146050	8698673445925	46	3
20117480050	8698673461611	20127148050	8698673445215	48	3
20117500050	8698673446618	20127500050	8698673429734	50	3




CORTEX SCREW CONICAL HEAD SELFT Ø 5.0 MM

Stainless Steel		Titanium		Length	Qty
TST Code	Barcode	TST Code	Barcode		
20110120050	8698673481497	20120120050	8698673481480	12	3
20110140050	8698673461796	20120140050	8698673475939	14	3
20110160050	8698673446700	20120160050	8698673429369	16	3
20110180050	8698673446717	20120180050	8698673429376	18	3
20110200050	8698673446724	20120200050	8698673429383	20	3
20110220050	8698673446731	20120220050	8698673429390	22	3
20110240050	8698673446748	20120240050	8698673429406	24	3
20110260050	8698673446755	20120260050	8698673429413	26	3
20110280050	8698673446762	20120280050	8698673429420	28	3
20110300050	8698673446779	20120300050	8698673429437	30	3
20110320050	8698673446786	20120320050	8698673429444	32	3
20110340050	8698673446793	20120340050	8698673429451	34	3
20110360050	8698673446809	20120360050	8698673429468	36	3
20110380050	8698673446816	20120380050	8698673429475	38	3
20110400050	8698673446823	20120400050	8698673429482	40	3
20110420050	8698673480643	20120420050	8698673480605	42	3
20110440050	8698673480650	20120440050	8698673480612	44	3
20110460050	8698673480667	20120460050	8698673480629	46	3
20110480050	8698673480674	20120480050	8698673480636	48	3
20110500050	8698673446847	20120500050	8698673429505	50	3



LOCK CORTICO CANCEL CONICAL HEAD SELFT Ø 5.5 MM

Stainless Steel		Titanium		Length	Qty
TST Code	Barcode	TST Code	Barcode		
20117500055	8698673446199	20127500055	8698673438194	50	3
20117550055	8698673446205	20127550055	8698673438200	55	3
20117600055	8698673446212	20127600055	8698673438217	60	3
20117650055	8698673446229	20127650055	8698673438224	65	3
20117700055	8698673446236	20127700055	8698673438231	70	3
20117750055	8698673446243	20127750055	8698673438248	75	3
20117800055	8698673446250	20127800055	8698673438255	80	3
20117850055	8698673446267	20127850055	8698673438262	85	3
20117900055	8698673446274	20127900055	8698673438279	90	3
20117950055	8698673446281	20127950055	8698673438286	95	3
20117100055	8698673446298	20127100055	8698673438149	100	3
20117105055	8698673446304	20127105055	8698673438156	105	3
20117110055	8698673446311	20127110055	8698673438163	110	3



CORTICO CANCELLOUS CONICAL HEAD SELFT Ø 5.5 MM

Stainless Steel		Titanium		Length	Qty
TST Code	Barcode	TST Code	Barcode		
20110500055	8698673446342	20120500055	8698673437883	50	2
20110550055	8698673446359	20120550055	8698673437890	55	2
20110600055	8698673446366	20120600055	8698673437906	60	2
20110650055	8698673446373	20120650055	8698673437913	65	2
20110700055	8698673446380	20120700055	8698673437920	70	2
20110750055	8698673446397	20120750055	8698673437937	75	2
20110800055	8698673446403	20120800055	8698673437944	80	2
20110850055	8698673446410	20120850055	8698673437951	85	2
20110900055	8698673446427	20120900055	8698673437968	90	2
20110950055	8698673446434	20120950055	8698673437975	95	2
20110100055	8698673446441	20120100055	8698673437982	100	2
20110105055	8698673446458	20120105055	8698673437999	105	2
20110110055	8698673446465	20120110055	8698673438002	110	2

LW-01 SS / LW-02 Ti Cable Plates Sets

Plates & Instruments

No	TST Code	Barcode	Description	Quantity
	00000804100	8699931026023	CABLE GRIP 1. DESIGN TRAY	1
	36810003001	8699931027730	PROX. FEMUR LOCK. HOOK PLATE S.S. SMALL	1
	36810002002	8699931027747	PROX. FEMUR LOCK. HOOK PLATE S.S. MEDIUM	1
	36810000130	8699931027754	PROX. FEMUR LOCK. HOOK PLATE S.S. 130 MM (3 HOLES)	1
	36815000180	8699931027761	PROX. FEMUR LOCK. HOOK PLATE S.S. R 180MM(5 HOLES)	1
	36815000255	8699931027778	PROX. FEMUR LOCK. HOOK PLATE S.S. R 255MM(8 HOLES)	1
LW-01 SS	36816000180	8699931027792	PROX. FEMUR LOCK. HOOK PLATE S.S. L 180MM(5 HOLES)	1
	36816000255	8699931027808	PROX. FEMUR LOCK. HOOK PLATE S.S. L 255MM(8 HOLES)	1
	36817000040	8699931028218	DIAPHYSEAL LOCK. CABLE PLATE S.S. 4 HOLES	1
	36817000060	8699931028225	DIAPHYSEAL LOCK. CABLE PLATE S.S. 6 HOLES	1
	36817000080	8699931028232	DIAPHYSEAL LOCK. CABLE PLATE S.S. 8 HOLES	1
	36817000100	8699931028249	DIAPHYSEAL LOCK. CABLE PLATE S.S. 10 HOLES	1
	26810000075	8699931027839	CABLE CLAMP SCREW	30
	36810000004	8699931027822	CABLE CLAMP	30
	30110000002	8680858420636	2.0 MM CABLE+SWAGE CRNI	5
	00000804100	8699931026023	CABLE GRIP 1. DESIGN TRAY	1
	36820003001	8699931027655	PROX. FEMUR LOCK. HOOK PLATE SMALL	1
	36820002002	8699931027662	PROX. FEMUR LOCK. HOOK PLATE MEDIUM	1
	36820000130	8699931027679	PROX. FEMUR LOCK. HOOK PLATE 3 HOLES	1
	36825000180	8699931027686	PROX. FEMUR LOCK. HOOK PLATE 5 HOLES	1
	36825000255	8699931027693	PROX. FEMUR LOCK. HOOK PLATE 8 HOLES	1
LW-02 Ti	36826000180	8699931027716	PROX. FEMUR LOCK. HOOK PLATE 5 HOLES	1
	36826000255	8699931027723	PROX. FEMUR LOCK. HOOK PLATE 8 HOLES	1
	36827000040	8699931028256	DIAPHYSEAL LOCK. CABLE PLATE 4 HOLES	1
	36827000060	8699931028263	DIAPHYSEAL LOCK. CABLE PLATE 6 HOLES	1
	36827000080	8699931028270	DIAPHYSEAL LOCK. CABLE PLATE 8 HOLES	1
	36827000100	8699931028287	DIAPHYSEAL LOCK. CABLE PLATE 10 HOLES	1
	26820000075	8699931028102	CABLE CLAMP SCREW	30
	36820000004	8699931028058	CABLE CLAMP	30
	30110000002	8680858420636	2.0 MM CABLE+SWAGE CRNI	5
1	00000804200	8699931026030	CABLE GRIP 2. DESIGN TRAY	1
2	03006100050	8698673496392	CABLE CUTTER	1
3	03006100020	8698673495418	FEMORAL CABLE PASSER	1
4	03006100030	8699931023626	LARGE CABLE PASSER	1
5	02010101002	8698673493308	SOFT SCREW DRIVER QUICK LARGE	1
6	02020101035	8698673493322	SCREW DRIVER QUICK Ø 3.5 MM	1
7	02231030042	8698673484795	BONE DRILL 4.2X60X300 MM	2
7	23410250120	8698673453128	KIRSCHNER WIRE TROCAR POINT 2X250 MM	4
7	04550129310	8680858424801	K-WIRE TUBE Ø12XØ9X310 MM	1
8	02005000100	8698673497047	DEPTH GAUGE - LOWER EXTREMITY PLATE	1
9	03501000042	8699931011005	CONIC THR.DRILL SLEEVE 4.2MM	1
10	08050000042	8699931026733	DRILL GUIDE Ø 4.2 MM	1
11	00000804300	8699931026047	CABLE GRIP 3. DESIGN TRAY	1
12	03006100040	8698673496934	CRIMP TOOL	1
13	03006100010	8698673496415	CABLE TENSIONER	1
14	00000804400	8699931028119	CABLE GRIP 4. DESIGN TRAY	1
15	08050000080	8699931027907	TROCHANTERIC GRIP POSITIONER	1
16	08060000005	8699931027914	TROCHANTERIC GRIP IMPACTOR	1
17	01193002009	8699931028126	BONE HAMMER - MEDIUM	1
	00000810100	8699931026054	CABLE GRIP SCREW TRAY	1
	00560270150	8699931010275	CONTAINER 560X270X150 MM	1

Instruments

