

UP-03 Small Head Distal Radius Plates Set

WRIST FUSION LOCKING PLATES*

Indications: This system is indicated for wrist arthrodesis and also for specific indications including posttraumatic arthritis of the joints of the wrist, rheumatoid wrist deformities requiring restoration, complex carpal instability, sequelae of septic arthritis of the wrist, severe unremitting wrist pain related to motion, brachial plexus nerve palsies, tumor resection and spastic deformities.

Plate Profile: Thickness: 2.5 mm, Width: 8.0 mm, Distance between holes: 12.0 mm.

TST Code	Barcode	Description	Material	Length (mm)
30320710001	8699931010152	WRIST FUSION LOCKING PLATE BEND STANDARD	Ti	142
30320710002	8699931010169	WRIST FUSION LOCKING PLATE BEND SHORT	Ti	130
30320700003	8699931010176	WRIST FUSION LOCKING PLATE STRAIGHT	Ti	112



SMALL HEAD DISTAL ULNA PLATES*

Indications: Fractures of the distal ulna which result in instability of the distal radioulnar joint and fractures of the ulna head where the articular surface is either displaced, rotated or tilted. Comminuted extra-articular fractures of the ulnar neck threatening stable congruency of the distal radioulnar joint.

Plate Profile: Thickness: 1.5 mm, Width: 6.0 mm.

TST Code	Barcode	Description	Material	Length (mm)
34828100005	8699931029635	SMALL HEAD DISTAL ULNA LOCK PLATE 5 HOLES	Ti	36
34828100006	8699931029642	SMALL HEAD DISTAL ULNA LOCK PLATE 6 HOLES	Ti	46
34828100007	8699931029659	SMALL HEAD DISTAL ULNA LOCK PLATE 7 HOLES	Ti	51



Recommended Screws 2.7 mm Small Head, 2.7 mm Small Head Lock, 3.5 mm Cortical, 3.5 mm Lock. Cortical