

UP-13
INTERCARPAL
FUSION
PLATES
SET

UP-13 Intercarpal Fusion Plates Set

INTERCARPAL FUSION PLATE Ti

The Mini Intercarpal Fusion Plate is an innovative implant for partial midcarpal arthrodesis of the wrist (Four-Corner-Fusion), for pseudo-arthrosis of the Os Scaphoid and aged SL-band ruptures (SLAC and SNAC wrist). This therapy has established itself as a treatment strategy. The application of it makes a compression between the merged carpal elements, hereby benefits a complete arthrodesis. The Sliding Plate offers the possibility to apply the compression centrally by means of turning the implant through sliding holes which are arranged in a circular way.

Indications: Posttraumatic mediocarpal arthritis following scaphoid non-union (SNAC) or scapholunate ligament rupture (SLAC) and mediocarpal instability (final stage).

Advantages:

Minimal outer diameter in three different sizes. Plate sizes; 15mm, 16mm and 17mm. Slight panel thickness/panel height 2.8mm.

Low profile circular plate with smooth surface and rounded edges to minimize soft tissue irritation.

Kirschner wire holes for temporary fixation and rotatory compression holes.

Fixable with 8 screws. Narrow screw diameter (2.4mm).

Biocompatible titanium alloy. Innovative instrumentation.

Variable screw angle (+/- 8°).

INTERCARPAL FUSION PLATE

TST Code	Barcode	Description	Size	Qty
36020000150	8699931020724	INTERCARPAL FUSION PLATE	Ø15 MM	2
36020000160	8699931020779	INTERCARPAL FUSION PLATE	Ø16 MM	2
36020000170	8699931020731	INTERCARPAL FUSION PLATE	Ø17 MM	2



SMALL HEAD SCREWS Ø 2.4 MM

TST Code	Barcode	Length	Qty
24620100024	8699931015218	10 MM	8
24620120024	8699931006728	12 MM	8
24620140024	8699931006735	14 MM	8
24620160024	8699931006742	16 MM	8
24620180024	8699931006759	18 MM	8



UP-13 Intercarpal Fusion Plates Set

Instruments

No	TST Code	Barcode	Description	Qty
1	00000403100	8699931024517	INTERCARPAL FUSION 1.DESIGN TRAY	1
2	00000757025	8699931024654	INTERCARPAL FUSION PLATE & SCREW BOX	1
3	02015100015	8699931021547	UNIVERSAL SCREW DRIVER HANDLE	1
4	02025010170	8698673485648	SCREW DRIVER HEXAGONAL Ø1.7	1
5	03650000015	8699931023855	INTERCARPAL FUSION PLATE REAMER 15 MM	1
6	03650000016	8699931023862	INTERCARPAL FUSION PLATE REAMER 16 MM	1
7	03650000017	8699931023879	INTERCARPAL FUSION PLATE REAMER 17 MM	1
8	03650000915	8699931023893	INTERCARPAL FUSION PLATE REAMER GUIDE 15 MM	1
9	03650000916	8699931023909	INTERCARPAL FUSION PLATE REAMER GUIDE 16 MM	1
10	03650000917	8699931023916	INTERCARPAL FUSION PLATE REAMER GUIDE 17 MM	1
11	03650000010	8699931023848	INTERCARPAL FUSION PLATE HOLDER	1
12	02350000020	8699931023602	BONE DRILL QUICK END Ø2 MM	2
13	02005000500	8698673497023	DEPTH GAUGE (0-50 MM)	1
14	00025502010	8698673448322	MOTOR QUICK ADAPTOR	1
15	04550129090	8699931032482	K-WIRE TUBE Ø12XØ9X90 MM	1
15	23410080115	8699931013429	KIRSCHNER WIRE TROCER POINT 1.5X80 MM	6
15	23410080112	8698673441231	KIRSCHNER WIRE TROCER POINT 1.2X80 MM	6
16	04010350070	8699931023923	L-ALLEN WRENCH MINI 3.5X70 MM	1
17	02005100020	8699931019094	DRILL GUIDE 2 MM (HANDLED)	1
18	03650001015	8699931028614	INTERCARPAL FUSION PLATE DRILL GUIDE 15 MM	1
19	03650001016	8699931028621	INTERCARPAL FUSION PLATE DRILL GUIDE 16 MM	1
20	03650001017	8699931028638	INTERCARPAL FUSION PLATE DRILL GUIDE 17 MM	1
	00270270080	8699931009019	CONTAINER 270X270X80 MM	1

