

tst | Orthopedic Implants

CS-03
HERBERT
SCREWS
SET

HERBERT SCREWS

Indications: Herbert Bone Screws feature short thread lengths and small diameters for optimum fixation of very small bone fragments. Manufactured from Titanium, the bone screws provide strength and biocompatibility during insertion. The double thread construct with a different pitch on distal and proximal ends, combines fracture reduction and compression to provide rigid internal fixation. Headless screws can be inserted through articular cartilage and buried below the bone surface to avoid interface with articulation or soft tissue. Leading threads penetrate bone faster than trailing threads, slowing penetration, reducing and compressing fragments during insertion.

TST Cannulated 2.5 mm Herbert Screws are indicated for fixation of fractures and nonunions of small bones and small bone arthrodeses, including scaphoid fractures; intra-articular fractures of the tarsals, metatarsals, carpals and metacarpals; bunionectomies and osteotomies; arthrodeses of small joints (e.g. phalanges); fractures of the patella, ulna and radial styloid.

3.0 mm Herbert Screws are intended for fixation of intra-articular and extra-articular fractures and nonunions of small bones and small bone fragments; arthrodeses of small joint; bunionectomies and osteotomies, including scaphoid and other carpal bones, metacarpals, tarsals, metatarsals, patella, ulnar styloid, capitellum, radial head and radial styloid.

4.3 mm Herbert Screws are indicated for fracture fixation, reconstruction, osteotomy and arthrodesis of various bones and bone fragments including joint fusions (arthrodeses) in the foot and fixation of intra-articular fractures of the humerus, femur and tibia.

TST CANNULATED HERBERT SCREW - Titanium

TST Code	Barcode	Description	Diameter (TIP X TAIL)	Lenght (TIP X TAIL)	Lenght TOTAL	Qty
20527001025	8699931028294	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	4 X 2 MM	10 MM	3
20527001225	8699931008746	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	4 X 2 MM	12 MM	3
20527001425	8699931008753	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	4 X 2 MM	14 MM	3
20527001625	8699931008760	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	5 X 2 MM	16 MM	3
20527001825	8699931008777	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	5 X 2 MM	18 MM	3
20527002025	8699931008814	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	5 X 2 MM	20 MM	3
20527002225	8699931008821	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	5 X 2 MM	22 MM	3
20527002425	8699931008838	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	6 X 2 MM	24 MM	3
20527002625	8699931008845	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	6 X 2 MM	26 MM	3
20527002825	8699931008852	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	6 X 2 MM	28 MM	3
20527003025	8699931008869	TST CANNULATED HERBERT SCREW	2.5 X 3.2 MM	6 X 2 MM	30 MM	3
20527001030	8698673469815	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	4 X 2.4 MM	10 MM	3
20527001230	8698673469822	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	4 X 2.4 MM	12 MM	3
20527001430	8698673469839	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	5 X 5 MM	14 MM	3
20527001630	8698673469846	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	5 X 5 MM	16 MM	3
20527001830	8698673469853	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	5 X 5 MM	18 MM	3
20527002030	8698673469860	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	5 X 5 MM	20 MM	3
20527002230	8698673469877	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	5 X 5 MM	22 MM	3
20527002430	8698673469884	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	5 X 5 MM	24 MM	3
20527002630	8698673469891	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	7 X 5 MM	26 MM	3
20527002830	8698673469907	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	7 X 5 MM	28 MM	3
20527003030	8698673469914	TST CANNULATED HERBERT SCREW	3.0 X 4.0 MM	7 X 5 MM	30 MM	3
20527001443	8698673482821	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	5 X 5 MM	14 MM	3
20527001643	8698673482838	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	5 X 5 MM	16 MM	3
20527001843	8698673482845	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	5 X 5 MM	18 MM	3
20527002043	8698673482852	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	5 X 5 MM	20 MM	3
20527002243	8698673482869	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	5 X 5 MM	22 MM	3
20527002443	8698673482876	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	5 X 5 MM	24 MM	3
20527002643	8698673482883	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	26 MM	3
20527002843	8698673482890	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	28 MM	3
20527003043	8698673482906	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	30 MM	3
20527003243	8698673482913	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	32 MM	3
20527003443	8698673482920	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	34 MM	3
20527003643	8698673482937	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	36 MM	3
20527003843	8698673482944	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	38 MM	3
20527004043	8698673482951	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	40 MM	3
20527004243	8698673482968	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	42 MM	3
20527004443	8698673482975	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	10 X 5 MM	44 MM	3
20527004643	8698673482982	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	46 MM	3
20527004843	8698673482999	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	48 MM	3
20527005043	8698673483002	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	50 MM	3
20527005243	8698673483019	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	52 MM	3
20527005443	8698673483026	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	54 MM	3
20527005643	8698673483033	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	56 MM	3
20527005843	8698673483040	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	58 MM	3
20527006043	8698673483057	TST CANNULATED HERBERT SCREW	4.3 X 5.3 MM	15 X 5 MM	60 MM	3



CS-03 Herbert Screws Set

Instruments

No	TST Code	Barcode	Description	Qty
1	0000901100	8699931017427	HERBERT SCREW 1. TRAY (260X260X65)	1
2	00040101400	8699931017434	HERBERT SCREW TRAY	1
3	02025000170	8698673469396	CANNULATED SCREW DRIVER QUICK HEX.TIP 1.7X80MM	1
4	02050100425	8699931005417	CANNULATED SCREW DRIVER QUICK TIP 2.5XØ5X130MM	1
5	02000452017	8699931029048	SCREW DRIVER QUICK TIP HEXAGONAL Ø1.7XØ4X52 MM	1
6	02020103025	8699931010930	SCREW DRIVER QUICK TIP 2.5XØ5X130 MM	1
7	02010101004	8699931021868	SOFT SCR. DRIVER CANNULATED QUICK SMALL	1
8	02065100002	8699931021899	Ø2 MM DRILL SLEEVE(HERBERT)	1
9	02065100025	8699931018905	DRILL SLEEVE GUIDE Ø 2.5 MM (HERBERT)	1
10	02065100032	8699931018073	Ø3.2 MM DRILL SLEEVE(HERBERT)	1
11	02065100045	8699931028898	Ø4.5 MM DRILL SLEEVE (HERBERT)	1
12	02065100001	8699931023565	Ø1 MM KIRSCHNER SLEEVE(HERBERT)	1
13	02065100015	8699931018059	Ø1.5 MM KIRSCHNER SLEEVE(HERBERT)	1
14	02065100055	8699931017946	Ø5.5 MM HERBERT SCREW SLEEVE	1
15	02050000130	8699931028928	LENGTH MEASURING DEVICE - CANN. HERBERT SCREW	1
16	03020000115	8698673493353	PENSET SMALL	1
17	04551000160	8699931028805	K-WIRE TUBE Ø10XØ8X160 MM	2
17	23410013012	8699931028904	KIRSCHNER WIRE TROCAR POINT 1X130 MM	5
17	23410013015	8699931028911	KIRSCHNER WIRE TROCAR POINT 1.5X130 MM	5
17	02062011201	8698673469402	CANNULATED DRILL 2X1.1X20X100 MM	1
17	02062511100	8699931028935	CANNULATED DRILL Ø 2.5X1.1X2.7X100 MM	1
17	02063216515	8698673485525	CANNULATED DRILL Ø 3.2X1.6X50X150	1
17	02064516150	8699931028959	CANNULATED DRILL Ø 4.5X1.6X15X150 MM	1
	00270270080	8699931009019	CONTAINER 270X270X80 MM	1

