


## UP-01 Upper Extremity Plates Set

## SCAPULA PLATES\*

**Indications:** Scapula body fractures, glenoid neck fractures, intra-articular glenoid fractures, scapula spine fractures and acromial fractures.

**Plate Profile:** Thickness: 3.0 mm, \*Glenoid: 2.5 mm, Width: 10.0 mm.

TST Code	Barcode	Description	Side	Material	Length (mm)
36526000101	8699931022155	SCAPULA LATERAL LOCKING LONG PLATE	L	Ti	100
36525000101	8699931022162	SCAPULA LATERAL LOCKING LONG PLATE	R	Ti	100
36526000201	8699931025149	SCAPULA LATERAL LOCKING SHORT PLATE	L	Ti	71
36525000201	8699931025132	SCAPULA LATERAL LOCKING SHORT PLATE	R	Ti	71
36526010200	8699931022216	SCAPULA MEDIAL LOCKING PLATE LONG	L	Ti	96 / 100
36525010200	8699931022223	SCAPULA MEDIAL LOCKING PLATE LONG	R	Ti	96 / 100
36526020200	8699931022230	SCAPULA MEDIAL LOCKING PLATE STANDARD	L	Ti	73 / 100
36525020200	8699931022247	SCAPULA MEDIAL LOCKING PLATE STANDARD	R	Ti	73 / 100
36526010020	8699931022292	SCAPULA SPINE LOCKING PLATE LONG	L	Ti	76
36525010020	8699931022308	SCAPULA SPINE LOCKING PLATE LONG	R	Ti	76
36526020020	8699931022315	SCAPULA SPINE LOCKING PLATE STANDARD	L	Ti	72
36525020020	8699931022322	SCAPULA SPINE LOCKING PLATE STANDARD	R	Ti	72
36526000010	8699931022117	GLENOID LOCKING PLATE	L	Ti	27
36525000010	8699931022124	GLENOID LOCKING PLATE	R	Ti	27



**Recommended Screws** 3.5 mm Cortical, 3.5 mm Lock. Cortical, 4.0 mm Cancellous, 4.0 mm Lock. Cancellous