

LW-22 Ankle Plates Set

ANATOMIC DISTAL FIBULA VAL PLATE

Indications: Anatomic Distal Fibula Plates are indicated for fractures, osteotomies and nonunions of the metaphyseal and diaphyseal region of the distal fibula.

Plate Profile: Thickness: 2.2 mm, Width: 11.0 mm, Distance between holes: 11.0 mm.

2.7 mm variable angled screws (VAL) allow up to 15° off axis angulation. VAL screws provide fixation options for a variety of fracture types.

TST Code	Barcode	Description	Side	Material
39325710008	8699931002478	ANATOMIC DIST.FIBULA LOCK. PLATE 15 HOLE	R	Ti
39325710010	8699931002485	ANATOMIC DIST.FIBULA LOCK. PLATE 17 HOLE	R	Ti
39325710012	8699931002492	ANATOMIC DIST.FIBULA LOCK. PLATE 19 HOLE	R	Ti
39326710008	8699931002515	ANATOMIC DIST.FIBULA LOCK. PLATE 15 HOLE	L	Ti
39326710010	8699931002522	ANATOMIC DIST.FIBULA LOCK. PLATE 17 HOLE	L	Ti
39326710012	8699931002539	ANATOMIC DIST.FIBULA LOCK. PLATE 19 HOLE	L	Ti

