


LW-22 Ankle Plates Set

CALCANEUS PLATES

Indications: Complicated calcaneal fractures, calcaneal osteotomies and severely comminuted fractures.

Plate Profile: Thickness 1.2 mm

TST Code	Barcode	Description	Material	Length (mm)
35125010000	8698673418622	LOWPROFILE CALCANEUS INTERLOCK. PLATE SHORT	Ti	58
35125020000	8698673418646	LOWPROFILE CALCANEUS INTERLOCK. PLATE STANDARD	Ti	67
35125030000	8698673418608	LOWPROFILE CALCANEUS INTERLOCK. PLATE LONG	Ti	75



Recommended Screw 3.5 mm Cortical, 3.5 mm Lock. Cortical, 4.0 mm Cancellous, 4.0 mm Lock. Cancellous