


## LW-22 Ankle Plates Set

# DISTAL FIBULA PLATES

**Indications:** Distal Fibula Locking Plates are intended for fixation of fractures, osteotomies and nonunions of the metaphyseal and diaphyseal region of the distal fibula.

**Plate Profile:** Thickness: 2.0 mm, Width: 10.5 mm, Distance between holes: 11.0 mm.

TST Code	Barcode	Description	Side	Material	Length (mm)
36525700008	8698673478022	LOW PROFILE DISTAL FIBULA LOCK.PLATE 8 HOLES	R	Ti	72
36525700010	8698673425804	LOW PROFILE DISTAL FIBULA LOCK.PLATE 10 HOLES	R	Ti	93
36525700012	8698673426771	LOW PROFILE DISTAL FIBULA LOCK.PLATE 12 HOLES	R	Ti	115
36526700008	8698673422438	LOW PROFILE DISTAL FIBULA LOCK.PLATE 8 HOLES	L	Ti	72
36526700010	8698673425927	LOW PROFILE DISTAL FIBULA LOCK.PLATE 10 HOLES	L	Ti	93
36526700012	8698673426764	LOW PROFILE DISTAL FIBULA LOCK.PLATE 12 HOLES	L	Ti	115



**Recommended Screw** 3.5 mm Cortical, 3.5 mm Lock. Cortical, 4.0 mm Cancellous, 4.0 mm Lock. Cancellous