

LW-36
LOWER
EXTREMITY
PLATES
SET


LW-36 Lower Extremity Plates Set

PROXIMAL FEMUR PLATES

Indications: Proximal Femur Plate is indicated for the treatment of fractures of the trochanteric region including simple intertrochanteric, reverse intertrochanteric, transverse trochanteric, complex multifragmentary and fractures with medial cortex instability. Proximal femur fractures with ipsilateral shaft fractures and metastatic proximal femur fractures. Proximal femur osteotomies, fractures in osteopenic bone, nonunions and malunions, basi/transcervical femoral neck fractures, subcapital femoral neck fractures and subtrochanteric femur fractures.

Plate Profile: Thickness: 6.5 mm, Width: 18.0 mm, Distance between holes: 21.5 mm

TST Code	Barcode	Description	Side	Material	Length (mm)
31625700007	8699931016178	PROX. FEMUR LOCKING PLATE CONICAL HOLES 7 H	R	Ti	155
31625700009	8699931016161	PROX. FEMUR LOCKING PLATE CONICAL HOLES 9 H	R	Ti	198
31625700011	8699931016192	PROX. FEMUR LOCKING PLATE CONICAL HOLES 11 H	R	Ti	241
31625700013	8699931016185	PROX. FEMUR LOCKING PLATE CONICAL HOLES 13 H	R	Ti	285
31626700007	8699931016215	PROX. FEMUR LOCKING PLATE CONICAL HOLES 7 H	L	Ti	155
31626700009	8699931016208	PROX. FEMUR LOCKING PLATE CONICAL HOLES 9 H	L	Ti	198
31626700011	8699931016239	PROX. FEMUR LOCKING PLATE CONICAL HOLES 11 H	L	Ti	241
31626700013	8699931016222	PROX. FEMUR LOCKING PLATE CONICAL HOLES 13 H	L	Ti	285



9 holes

Recommended Screws 5.0 mm Cortical, 5.0 mm Lock. Cortical, 5.5 mm Cancellous, 5.5 mm Lock. Cancellous

LW-36 Lower Extremity Plates Set

FEMUR PLATES

Indications: Osteosynthesis of femoral shaft fractures allowing buttressing, neutralization, tension banding, bridging, compression.

Plate Profile: Thickness: 6.5 mm, Width: 18.0 mm, Distance between holes: 21.0 mm

Cat. Code	Barcode	Description	Material	Length (mm)
30227100006	8698673438101	LC DCP CONICAL BOW LOCK FEMORAL PLATE 6 H	Ti	142
30227100007	8698673438118	LC DCP CONICAL BOW LOCK FEMORAL PLATE 7 H	Ti	168
30227100008	8698673438125	LC DCP CONICAL BOW LOCK FEMORAL PLATE 8 H	Ti	190
30227100009	8698673438132	LC DCP CONICAL BOW LOCK FEMORAL PLATE 9 H	Ti	212
30227100010	8698673438057	LC DCP CONICAL BOW LOCK FEMORAL PLATE 10 H	Ti	233
30227100011	8698673438064	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 11 H	Ti	254
30227100012	8698673438071	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 12 H	Ti	276
30227100013	8698673438088	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 13 H	Ti	298
30227100014	8698673438095	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 14 H	Ti	320



12 holes

LW-36 Lower Extremity Plates Set

MISS Distal Femur Plates

Indications: Distal shaft fractures, Periprosthetic fractures, supracondylar fractures and intra-articular fractures.

Plate Profile: Thickness: 6.0 mm, Width: 16.0 mm, Distance between holes: 20.5 mm.

TST Code	Barcode	Description	Side	Material	Length (mm)
35125710111	8698673428638	MISS LC DISTAL FEMUR PLATE 11 HOLES	R	Ti	152
35125710113	8698673428645	MISS LC DISTAL FEMUR PLATE 13 HOLES	R	Ti	193
35125710115	8698673428737	MISS LC DISTAL FEMUR PLATE 15 HOLES	R	Ti	234
35125710117	8698673428744	MISS LC DISTAL FEMUR PLATE 17 HOLES	R	Ti	275
35126710111	8698673428591	MISS LC DISTAL FEMUR PLATE 11 HOLES	L	Ti	152
35126710113	8698673428607	MISS LC DISTAL FEMUR PLATE 13 HOLES	L	Ti	193
35126710115	8698673428614	MISS LC DISTAL FEMUR PLATE 15 HOLES	L	Ti	234
35126710117	8698673428621	MISS LC DISTAL FEMUR PLATE 17 HOLES	L	Ti	275



17 holes

Recommended Screws 5.0 mm Cortical, 5.0 mm Lock. Cortical, 5.5 mm Cancellous, 5.5 mm Lock. Cancellous

LW-36 Lower Extremity Plates Set

MISS PROXIMAL TIBIA MEDIAL PLATES

Indications: The plate is indicated for metaphyseal fractures of the medial tibial plateau and split-type fractures of the medial tibial plateau.

Plate Profile: Thickness: 4.0 mm, Width: 15.5 mm, Distance between holes: 15.0 mm.

TST Code	Barcode	Description	Side	Material	Length (mm)
35025710207	8698673446007	MISS LC PROX. TIBIA MEDIAL PLATE 7 HOLES	R	Ti	89
35025710209	8698673446014	MISS LC PROX. TIBIA MEDIAL PLATE 9 HOLES	R	Ti	117
35025710211	8698673446021	MISS LC PROX. TIBIA MEDIAL PLATE 11 HOLES	R	Ti	145
35025710213	8698673446038	MISS LC PROX. TIBIA MEDIAL PLATE 13 HOLES	R	Ti	185
35025710215	8698673440227	MISS LC PROX. TIBIA MEDIAL PLATE 15 HOLES	R	Ti	225
35026710207	8698673444973	MISS LC PROX. TIBIA MEDIAL PLATE 7 HOLES	L	Ti	89
35026710209	8698673446045	MISS LC PROX. TIBIA MEDIAL PLATE 9 HOLES	L	Ti	117
35026710211	8698673446052	MISS LC PROX. TIBIA MEDIAL PLATE 11 HOLES	L	Ti	145
35026710213	8698673446069	MISS LC PROX. TIBIA MEDIAL PLATE 13 HOLES	L	Ti	185
35026710215	8698673440234	MISS LC PROX. TIBIA MEDIAL PLATE 15 HOLES	L	Ti	225



Recommended Screws 5.0 mm Cortical, 5.0 mm Lock. Cortical, 5.5 mm Cancellous, 5.5 mm Lock. Cancellous

LW-36 Lower Extremity Plates Set

MISS PROXIMAL TIBIA LATERAL PLATES

Indications: Split-type fractures of the lateral tibial plateau and lateral split fractures with associated depressions.

Plate Profile: Thickness: 4.0 mm, Width: 16.0 mm, Distance between holes: 14.0 mm..

TST Code	Barcode	Description	Side	Material	Length (mm)
35025710107	8698673444836	MISS LC PROXIMAL TIBIA LATERAL PLATE 7 HOLES	R	Ti	103
35025710109	8698673444843	MISS LC PROXIMAL TIBIA LATERAL PLATE 9 HOLES	R	Ti	133
35025710111	8698673438521	MISS LC PROXIMAL TIBIA LATERAL PLATE 11 HOLES	R	Ti	174
35025710113	8698673444850	MISS LC PROXIMAL TIBIA LATERAL PLATE 13 HOLES	R	Ti	215
35026710107	8698673444881	MISS LC PROXIMAL TIBIA LATERAL PLATE 7 HOLES	L	Ti	103
35026710109	8698673444898	MISS LC PROXIMAL TIBIA LATERAL PLATE 9 HOLES	L	Ti	133
35026710111	8698673444911	MISS LC PROXIMAL TIBIA LATERAL PLATE 11 HOLES	L	Ti	174
35026710113	8698673444935	MISS LC PROXIMAL TIBIA LATERAL PLATE 13 HOLES	L	Ti	215



13 holes

Recommended Screws 5.0 mm Cortical, 5.0 mm Lock. Cortical, 5.5 mm Cancellous, 5.5 mm Lock. Cancellous

LW-36 Lower Extremity Plates Set

TIBIA PLATES

Indications: Osteosynthesis of bones, such as femur or tibia, allowing buttressing, neutralization, tension banding, bridging and compression.

Plate Profile: Thickness: 5.0 mm, Width: 14.0 mm, Distance between holes: 21.0 mm

TST Code	Barcode	Description	Material	Length (mm)
30127000106	8699931013788	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 6 HOLES	Ti	138
30127000107	8699931013795	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 7 HOLES	Ti	159
30127000108	8699931013801	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 8 HOLES	Ti	180
30127000109	8699931013818	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 9 HOLES	Ti	201
30127000110	8699931013825	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 10 HOLES	Ti	222
30127000111	8699931013832	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 11 HOLES	Ti	243
30127000112	8699931013849	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 12 HOLES	Ti	264
30127000113	8699931013856	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 13 HOLES	Ti	285
30127000114	8699931013863	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 14 HOLES	Ti	306



8 holes

LW-36 Lower Extremity Plates Set

DISTAL FEMUR MEDIAL OSTEOTOMY PLATES*

Indications: Unicompartmental lateral gonarthrosis with valgus malalignment of the distal femur and idiopathic or posttraumatic valgus deformity of the distal femur. Distal femur medial plates are for closed-wedge osteotomies and have a fixed angle construct for stable fixation. The plates are available in right or left version and are designed according to the principles of the locking compression plate. In the distal section there are four threaded holes and the directions of them are adapted to the anatomy of the supracondylar femur.

Plate Profile: Thickness: 4.5 mm, Width: 15.0 mm, Distance between holes: 18.0 mm.

TST Code	Barcode	Description	Side	Material	Length (mm)
35125730008	8698673463295	LCP DIST.FEM.MED.OSTEOT.LOCK PLATE 8 HOLES	R	Ti	115
35126730008	8698673463318	LCP DIST.FEM.MED.OSTEOT.LOCK PLATE 8 HOLES	L	Ti	115



8 holes

LW-36 Lower Extremity Plates Set

Instruments



LOCKING CORTEX SCREW CONICAL HEAD SELFT Ti Ø 5 MM

TST Code	Barcode	Length	Qty
20127180050	8698673429604	18 mm	5
20127200050	8698673429611	20 mm	5
20127220050	8698673429628	22 mm	5
20127240050	8698673429635	24 mm	5
20127260050	8698673429642	26 mm	5
20127280050	8698673429659	28 mm	5
20127300050	8698673429666	30 mm	5
20127320050	8698673429673	32 mm	5
20127340050	8698673429680	34 mm	5
20127360050	8698673429697	36 mm	5
20127380050	8698673429703	38 mm	5
20127400050	8698673429710	40 mm	5
20127450050	8698673429727	45 mm	5
20127500050	8698673429734	50 mm	5
20127550050	8698673429741	55 mm	5
20127600050	8698673429758	60 mm	5
20127650050	8698673429765	65 mm	5
20127700050	8698673429772	70 mm	5



CORTEX SCREW CONICAL HEAD Ti Ø 5 MM

TST Code	Barcode	Length	Qty
20120180050	8698673429376	18 mm	4
20120200050	8698673429383	20 mm	4
20120220050	8698673429390	22 mm	4
20120240050	8698673429406	24 mm	4
20120260050	8698673429413	26 mm	4
20120280050	8698673429420	28 mm	4
20120300050	8698673429437	30 mm	4
20120320050	8698673429444	32 mm	4
20120340050	8698673429451	34 mm	4
20120360050	8698673429468	36 mm	4
20120380050	8698673429475	38 mm	4
20120400050	8698673429482	40 mm	4
20120450050	8698673429499	45 mm	4
20120500050	8698673429505	50 mm	4
20120550050	8698673429512	55 mm	4
20120600050	8698673429529	60 mm	4
20120650050	8698673429536	65 mm	4
20120700050	8698673429543	70 mm	4



LOCK CORTICO CANCEL. CONICAL HEAD SELFT Ti Ø 5.5 MM

TST Code	Barcode	Length	Qty
20127500055	8698673438194	50 mm	4
20127550055	8698673438200	55 mm	4
20127600055	8698673438217	60 mm	4
20127650055	8698673438224	65 mm	4
20127700055	8698673438231	70 mm	4
20127750055	8698673438248	75 mm	4
20127800055	8698673438255	80 mm	4
20127850055	8698673438262	85 mm	4
20127900055	8698673438279	90 mm	4
20127950055	8698673438286	95 mm	4
20127100055	8698673438149	100 mm	4
20127105055	8698673438156	105 mm	4



CORTICO CANCELLOUS CONICAL HEAD SELFT Ti Ø 5.5 MM

TST Code	Barcode	Length	Qty
20120500055	8698673437883	50 mm	4
20120550055	8698673437890	55 mm	4
20120600055	8698673437906	60 mm	4
20120650055	8698673437913	65 mm	4
20120700055	8698673437920	70 mm	4
20120750055	8698673437937	75 mm	4
20120800055	8698673437944	80 mm	4
20120850055	8698673437951	85 mm	4
20120900055	8698673437968	90 mm	4
20120950055	8698673437975	95 mm	4
20120100055	8698673437982	100 mm	4
20120105055	8698673437999	105 mm	4



CANNULATED LOCK.CANCEL.SCREW Ti Ø 5.5 MM

TST Code	Barcode	Length	Qty
20327950055	8680858405121	50 mm	2
20327955055	8680858405138	55 mm	2
20327960055	8680858405145	60 mm	2
20327965055	8680858405152	65 mm	2
20327970055	8680858405169	70 mm	2
20327975055	8680858405176	75 mm	2
20327980055	8680858405183	80 mm	2
20327985055	8680858405190	85 mm	2
20327990055	8680858405206	90 mm	2
20327995055	8680858405213	95 mm	2
20327910055	8680858405220	100 mm	2
20327910555	8680858405237	105 mm	2



CANNULATED PARTIAL THREAD SCREW Ti Ø 5.5 MM

TST Code	Barcode	Length	Qty
20328950055	8680858405244	50 mm	2
20328955055	8680858405251	55 mm	2
20328960055	8680858405268	60 mm	2
20328965055	8680858405275	65 mm	2
20328970055	8680858405282	70 mm	2
20328975055	8680858405299	75 mm	2
20328980055	8680858405305	80 mm	2
20328985055	8680858405312	85 mm	2
20328990055	8680858405329	90 mm	2
20328995055	8680858405336	95 mm	2
20328910055	8680858405343	100 mm	2
20328910555	8680858405350	105 mm	2



SCREW WASHER Ti Ø 13 MM

TST Code	Barcode	Diameter	Qty
20820000013	8698673428683	13	5

LW-36 Lower Extremity Plates Set

Plates & Instruments

No	TST Code	Barcode	Description	Qty
1	00000201310	8680858443925	LOWER EXTREMITY PLATES CONICAL 1. DESIGN TRAY	1
	31625700007	8699931016178	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI R 7 H	1
	31625700009	8699931016161	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI R 9 H	1
	31625700011	8699931016192	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI R 11 H	1
	31625700013	8699931016185	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI R 13 H	1
	31626700007	8699931016215	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI L 7 H	1
	31626700009	8699931016208	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI L 9 H	1
	31626700011	8699931016239	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI L 11 H	1
	31626700013	8699931016222	PROX. FEMUR LOCKING PLATE CONICAL HOLES TI L 13 H	1
2	03160000001	8698673462250	MISS PROXIMAL FEMUR GUIDE R	1
3	03160000002	8698673462243	MISS PROXIMAL FEMUR GUIDE L	1
4	03160000007	8698673463783	MISS PROXIMAL FEMUR GUIDE IMPACTOR	1
5	03160000905	8699931024814	MISS PROXIMAL FEMUR GUIDE LOCKING CONICAL SCREW	1
6	08360045235	8699931030792	MISS PROX.FEMUR REDUC.DEVICE Ø 4.5X235MM	1
7	04050000008	8699931028379	WRENCH Ø 8 MM	1
8	00000201320	8680858443932	LOWER EXTREMITY PLATES CONICAL 2. DESIGN TRAY	1
9	02000415035	8699931029260	SCREW DRIVER QUICK HEX. TIP Ø 3.5X150 MM	1
10	02264260250	8699931030822	BONE DRILL 4.2X60X250MM CALIBR.FOR MISS PROX.FEMUR	2
11	02264560250	8699931030839	BONE DRILL 4.5X60X250MM CALIBR.FOR MISS PROX.FEMUR	1
12	02064580300	8699931032420	CANNULATED DRILL BIT Ø 4.5X2.1X45X250 MM	1
13	03501000042	8699931011005	CONIC THR.DRILL SLEEVE 4.2MM	2
14	03501000045	8699931030778	CONIC THR.DRILL SLEEVE 4.5MM	1
15	02006000100	8680858410972	DEPTH GAUGE - MISS LOWER EXT. OVER PLATE	1
16	01191000040	8699931008739	SCREW DRIVER TORQUE 3NM	1
17	02010101002	8698673493308	SOFT SCREW DRIVER QUICK LARGE	1
18	03160000004	8698673462267	MISS PROXIMAL FEMUR KIRSCHNER GUIDE	1
19	02055000335	8698673495715	CANNULATED SCREWDRIVER Ø 3.5 MM	1
20	03150020200	8699931030846	MEASURING DEVICE FOR K-WIRE 2X200 MM (PROX. FEMUR)	1
21	04551080210	8699931030808	K-WIRE TUBE Ø10XØ8X210 MM	1
21	23412200020	8698673452718	KIRSCHNER WIRE TREADED POINT 2X200 MM	2
21	23410200120	8698673453111	KIRSCHNER WIRE TROCER POINT 2X200 MM	2
	00000201330	8680858443949	LOWER EXTREMITY PLATES CONICAL 3. DESIGN TRAY	1
	30227100006	8698673438101	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 6 HOLES	1
	30227100007	8698673438118	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 7 HOLES	1
	30227100008	8698673438125	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 8 HOLES	1
	30227100009	8698673438132	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 9 HOLES	1
	30227100010	8698673438057	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 10 HOLES	1
	30227100011	8698673438064	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 11 HOLES	1
	30227100012	8698673438071	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 12 HOLES	1
	30227100013	8698673438088	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 13 HOLES	1
	30227100014	8698673438095	LC DCP CONICAL BOW LOCK FEMORAL PLATE TI 14 HOLES	1
	35125710111	8698673428638	MISS LC DISTAL FEMUR PLATE TI R 11 HOLES	1
	35125710113	8698673428645	MISS LC DISTAL FEMUR PLATE TI R 13 HOLES	1
	35125710115	8698673428737	MISS LC DISTAL FEMUR PLATE TI R 15 HOLES	1
	35125710117	8698673428744	MISS LC DISTAL FEMUR PLATE TI R 17 HOLES	1
	35126710111	8698673428591	MISS LC DISTAL FEMUR PLATE TI L 11 HOLES	1
	35126710113	8698673428607	MISS LC DISTAL FEMUR PLATE TI L 13 HOLES	1
	35126710115	8698673428614	MISS LC DISTAL FEMUR PLATE TI L 15 HOLES	1
	35126710117	8698673428621	MISS LC DISTAL FEMUR PLATE TI L 17 HOLES	1
	00000201340	8680858443956	LOWER EXTREMITY PLATES CONICAL 4. DESIGN TRAY	1
	35025710107	8698673444836	MISS LC PROXIMAL TIBIA LATERAL PLATE R 7 HOLES	1
	35025710109	8698673444843	MISS LC PROXIMAL TIBIA LATERAL PLATE R 9 HOLES	1
	35025710111	8698673438521	MISS LC PROXIMAL TIBIA LATERAL PLATE R 11 HOLES	1
	35025710113	8698673444850	MISS LC PROXIMAL TIBIA LATERAL PLATE R 13 HOLES	1
	35026710107	8698673444881	MISS LC PROXIMAL TIBIA LATERAL PLATE L 7 HOLES	1
	35026710109	8698673444898	MISS LC PROXIMAL TIBIA LATERAL PLATE L 9 HOLES	1
	35026710111	8698673444911	MISS LC PROXIMAL TIBIA LATERAL PLATE L 11 HOLES	1
	35026710113	8698673444935	MISS LC PROXIMAL TIBIA LATERAL PLATE L 13 HOLES	1
	35025710207	8698673446007	MISS LC PROXIMAL TIBIA MEDIAL PLATE R 7 HOLES	1
	35025710209	8698673446014	MISS LC PROXIMAL TIBIA MEDIAL PLATE R 9 HOLES	1
	35025710211	8698673446021	MISS LC PROXIMAL TIBIA MEDIAL PLATE R 11 HOLES	1
	35025710213	8698673446038	MISS LC PROXIMAL TIBIA MEDIAL PLATE R 13 HOLES	1
	35025710215	8698673440227	MISS LC PROXIMAL TIBIA MEDIAL PLATE R 15 HOLES	1
	35026710207	8698673444973	MISS LC PROXIMAL TIBIA MEDIAL PLATE L 7 HOLES	1
	35026710209	8698673446045	MISS LC PROXIMAL TIBIA MEDIAL PLATE L 9 HOLES	1
	35026710211	8698673446052	MISS LC PROXIMAL TIBIA MEDIAL PLATE L 11 HOLES	1
	35026710213	8698673446069	MISS LC PROXIMAL TIBIA MEDIAL PLATE L 13 HOLES	1
	35026710215	8698673440234	MISS LC PROXIMAL TIBIA MEDIAL PLATE L 15 HOLES	1
	30127000106	8699931013788	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 6 HOLES	1
	30127000107	8699931013795	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 7 HOLES	1
	30127000108	8699931013801	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 8 HOLES	1
	30127000109	8699931013818	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 9 HOLES	1
	30127000110	8699931013825	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 10 HOLES	1
	30127000111	8699931013832	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 11 HOLES	1
	30127000112	8699931013849	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 12 HOLES	1
	30127000113	8699931013856	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 13 HOLES	1
	30127000114	8699931013863	LC DCP INTERLOCKING CONICAL TIBIAL PLATES 14 HOLES	1
	00000205500	8699931023510	MISS FEMUR-TIBIA SCREW DESIGN TRAY	1
	00560270170	8699931010787	CONTAINER 560X270X170 MM	1

LW-36 Lower Extremity Plates Set

Instruments

