

tst | Orthopedic Implants

LW-28
OWO
PLATES
SET


LW-28 OWO Plates Set

OWO PLATES (HIGH TIBIAL OPEN WEDGE OSTEOTOMY PLATES)

Indications: High Tibial Open Wedge Osteotomy is a commonly used surgery method especially for young and active patients who have Medial Osteoarthritis with pain and Varus Deformity on their knees. The aim is the acute correction of the mechanical axis from the injured medial plateau towards the lateral plateau.

Plate Profile: Thickness: 3.0 mm.

TST Code	Barcode	Description	Side	Material	Qty
34325000205	8698673419315	OWO ANATOMIC PLATE	R	Ti	4
34326000205	8698673419322	OWO ANATOMIC PLATE	L	Ti	4
34420000800	8699931029802	OWO BLOCK 8 MM	-	Ti	2
34420001000	8699931029819	OWO BLOCK 10 MM	-	Ti	2
34420001250	8699931029826	OWO BLOCK 12.5 MM	-	Ti	2
34420001500	8699931029833	OWO BLOCK 15 MM	-	Ti	2
21620080050	8698673419391	OWO BLOCK SCREW	-	Ti	16



LOCKING CORTICAL SCREW Ti Ø 4.5 MM

TST Code	Barcode	Length	Qty
20227280045	8698673417854	28	2
20227300045	8698673417861	30	2
20227320045	8698673417878	32	2
20227340045	8698673417885	34	2
20227360045	8698673417892	36	2
20227380045	8698673417908	38	2
20227400045	8698673417915	40	2
20227420045	8698673417922	42	2
20227440045	8698673417939	44	2
20227460045	8698673417946	46	2
20227480045	8698673417953	48	2
20227500045	8698673417960	50	2
20227520045	8698673417977	52	2



CORTICAL SCREW Ti Ø 4.5 MM

TST Code	Barcode	Length	Qty
20220280045	8698673417359	28	2
20220300045	8698673417366	30	2
20220320045	8698673417373	32	2
20220340045	8698673417380	34	2
20220360045	8698673417397	36	2
20220380045	8698673417403	38	2
20220400045	8698673417410	40	2
20220420045	8698673417427	42	2
20220440045	8698673417434	44	2
20220460045	8698673417441	46	2
20220480045	8698673417458	48	2
20220500045	8698673417465	50	2
20220520045	8698673417472	52	2



LOCKING CANCELLOUS SCREW Ti Ø 6.5 MM

TST Code	Barcode	Length	Qty
20427300065	8698673416956	30	3
20427350065	8698673416963	35	3
20427400065	8698673416970	40	3
20427450065	8698673416987	45	3
20427500065	8698673416994	50	3
20427550065	8698673417007	55	3
20427600065	8698673417014	60	3
20427650065	8698673417021	65	3
20427700065	8698673417038	70	3
20427750065	8698673417045	75	3
20427800065	8698673417052	80	3



CANCELLOUS SCREW Ti Ø 6.5 MM

TST Code	Barcode	Length	Qty
20420300065	8698673421790	30	2
20420350065	8698673421806	35	2
20420400065	8698673421813	40	2
20420450065	8698673421820	45	2
20420500065	8698673421837	50	2
20420550065	8698673421844	55	2
20420600065	8698673421851	60	2
20420650065	8698673421868	65	2
20420700065	8698673421875	70	2
20420750065	8698673421882	75	2
20420800065	8698673421899	80	2

LW-28 OWO Plates Set

Instruments

No	TST Code	Barcode	Description	Quantity
1	0000805100	8699931027969	OWO PLATE 1.DESIGN TRAY	1
	34325000205	8698673419315	OWO ANATOMIC PLATES (Right, Left)	8
	34420000800	8699931029802	OWO BLOCK (8, 10, 12.5, 15 MM)	8
	21620080050	8698673419391	OWO BLOCK SCREW	16
2	08060000241	8698673494817	PLATE GUIDE R FOR OWO	1
3	08060000244	8698673494824	PLATE GUIDE L FOR OWO	1
4	02005100035	8699931026689	THREADED DRILL GUIDE Ø3.5 MM	1
5	02005100045	8699931026696	THREADED DRILL GUIDE Ø4.5 MM	1
6	03436025000	8698673494930	OSTEOTOME 25 MM	1
7	03436035000	8698673494909	OSTEOTOME 35 MM	1
8	03310016060	8680858409006	CHISEL BLADE 16X60 MM STRAIGHT	1
8	03310011060	8680858409013	CHISEL BLADE 11X60 MM STRAIGHT	1
9	03381000001	8699931023770	CHISEL HANDLE	1
10	04551080210	8699931030808	K-WIRE TUBE Ø10XØ8X210 MM	3
10	22310160035	8698673458826	BONE DRILL 3.5X60X160 MM	2
10	22310200045	8698673450493	BONE DRILL 4.5 X 200 MM	2
10	23410200115	8698673452992	KIRSCHNER WIRE TROCER POINT 1.5X200 MM	3
10	23410200120	8698673453111	KIRSCHNER WIRE TROCER POINT 2X200 MM	10
10	23410201025	8680858419524	KIRSCHNER WIRE TROCER POINT 2.5X200 MM	5
11	02005000100	8698673497047	DEPTH GAUGE - LOWER EXTREMITY PLATE	1
12	0000805200	8699931027952	OWO PLATE 2.DESIGN TRAY	1
13	08050000145	8698673497108	DISTRACTOR OSTEOTOME	1
14	02020103025	8699931010930	SCREW DRIVER QUICK TIP 2.5XØ5X130 MM	1
15	02010101002	8698673493308	SOFT SCREW DRIVER QUICK LARGE	1
16	02000415035	8699931029260	SCREW DRIVER QUICK HEX. TIP Ø 3.5X150 MM	1
17	03435000080	8698673494893	AUTOMATIC DISTRACTOR	1
18	01193001009	8699931028140	BONE HAMMER - LARGE	1
19	03435000030	8698673493513	WEDGE DOUBLE OSTEOTOME	1
20	08040000151	8698673494800	PLATE IMPACTOR	1
	0000805400	8699931027938	OWO PLATE SCREW TRAY	1
	00560270125	8699931010923	CONTAINER 560X270X125 MM	1

