


MULTIAXIAL DISTAL RADIUS WRIST FIXATOR



Another type of fixator for the hand-wrist application is the “Multiaxial Distal Radius Hand-Wrist fixator”.

It is composed of a central located C shaped body and bars which are connected with the ball-joint structure. Pin clamps can move on the bars for compression or distraction.

The C body of the fixator allows to view the fracture area.

CODE	BARCODE	DESCRIPTION	TYPE	LENGTH (mm)
40640020013	8698673456730	MULTIAXIAL DISTAL RADIUS WRIST FIXATOR	C	230

