

PelviPEDAL  
HIP SPICA TABLE

**PelviPEDAL**  
hip spica cast table


# PelviPEDAL HIP SPICA TABLE

PelviPedal Hip Spica Cast Table is designed and manufactured according to the specifications of the pediatric hip surgery. It is adjustable, efficient, user friendly, and provides both stability and maneuverability while applying a hip spica cast in the treatment of pediatric application. PelviPedal Hip Spica Cast Table is specifically tailored for child-care specialists working within a professional medical setting. Whether in clinics, children's hospitals, emergency rooms, general hospitals, the use of PelviPedal Hip Spica Cast Table will ease the treatment of pediatric fractures, osteotomy, and dislocations. It provides a stable and adjustable platform that is easily stowed and cleaned after use. Constructed of stainless steel with radiolucent platform, Simple and adjustable height, and length, Ideal for full body casts and hip spicas. PelviPedal also has Child's Conversion Kit.

Cat. Code	Barcode	DESCRIPTION
01040000000	8699931013177	PELVIPEDAL HIP SPICA CAST TABLE



# PeLviPEDAL HIP SPICA TABLE



PeLviPEDAL Hip Spica Table

